

ARJOHUNTLEIGH

GETINGE GROUP

AUTOLOGIC 200 MATTRESS REPLACEMENT



The *Auto logic*® 200 mattress replacement is part of the new logic range from ArjoHuntleigh.

The *Auto logic* 200 is recommended for round-the-clock prevention and management of up to grade 4 pressure ulcers with frequent monitoring, nursing intervention and repositioning.

Forget about adjusting pump settings during installation. For optimal pressure

area management of very high risk patients, the all-new *Auto logic* pump uses Self Set Technology (SST) to control the mattress pressure. In both Active (Alternating) and Reactive (Constant Lower Pressure) therapy modes, an advanced microprocessor system regularly assesses the body mass distribution of patients and readjusts cell air pressures to suit their individual requirements.

...with people in mind

www.ArjoHuntleigh.com

KEY FEATURES



- Rapid deflation to assist resuscitation
- Simple one handed operation
- Positioned at head end



Provides a temporary stable support surface to assist patient transfer and other nursing procedures




- Centralises patient on the mattress
- Increases stability of edge to aid transfer off and on the surface



Providing leading pressure management therapy options:


- Active (Alternating) mode
- Reactive (Constant Lower Pressure) mode




Allows continuity of chosen therapy for a minimum of 8 hours during patient transport or in the event of a power cut



- Simple operation via disconnection or tube-set from the pump
- Provides patient support for up to 12 hours



- Helps reduce humidity at the patient - mattress interface



- Standard *Dartex*® Cover
- *Advantex*® Cover - lower friction properties
- *eVENT*® Fabric Cover - lower friction and moisture vapour permeability

ArjoHuntleigh is a branch of Arjo Ltd. Med. AB. Only ArjoHuntleigh designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by ArjoHuntleigh. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice. © and ™ are trademarks belonging to the ArjoHuntleigh group of companies. © ArjoHuntleigh, 2011.

ArjoHuntleigh

310-312, Dallow Road Luton Bedfordshire LU1 1TD United Kingdom • Phone: +44 1582 745 800 • Fax: +44 1582 745 866

E-mail: international@arjohuntleigh.com • www.ArjoHuntleigh.com